



Niche NHS
Our Focus Is You

DAY CARE CHAIRS





Mobile, reclining, battery powered day chair

Ref. Stirling



Fabric choice

- Modern, stylish day chair with comfortably padded upholstery
- Battery powered recline to the most comfortable, practical position by using the recessed button outside the right arm (when seated) – charged overnight it gives approx. 60/70 full cycles per day
- Wipe clean upholstery on all areas which come into contact with the occupant - smooth surfaces without buttons or piping
- Fitted with four swivel, total lock brake castors for safety - no wires/ trip hazards for staff or users. These must be braked when the user is sitting down/alighting or if the chair has been left unoccupied
- Full width pushing handle fitted at the rear to assist when moving the unoccupied chair

Please note: These are heavy day chairs and are not for portering. They must only be used indoors on level and dry surfaces. They must not be moved when occupied.



Fully reclined



Semi recline



Easy to use button on outside of right arm when seated

Height	1060mm (42")
Width	810mm (32")
Depth	970mm (38")
Seat Height	490mm (19")
Seat Width	540mm (21")
Seat Depth	490mm (19")
Arm Height	200mm (8")
Back Height	730mm (29")
Max User Weight	175kg (25st)
Product Weight	50kg



Mobile, reclining, battery powered day chair

Ref. Hanley



Fabric choice

- Gives freedom to move the unoccupied chair as you like
- Padded upholstery with spoon shape backrest for napping
- Battery powered recline to the most comfortable, practical position by using the recessed button outside the right arm (when seated) - charged overnight it gives approx. 60/70 full cycles per day
- Fitted with four swivel, total lock brake castors for safety - no wires/trip hazards for staff or users. These must be braked when the user is sitting down/alighting or if the chair has been left unoccupied
- Full width pushing handle fitted at the rear to assist when moving the unoccupied chair

Please note: These are heavy day chairs and are not for portering. They must only be used indoors on level and dry surfaces. They must not be moved when occupied.



Wide pushing handle



Easy to use button



Height	1110mm (44")
Width	740mm (29")
Depth (Seat Retracted)	925mm (36")
Seat Height	540mm (21")
Seat Width	549mm (21")
Seat Depth	490mm (19")
Arm Height	180mm (7")
Back Height	740mm (29")
Max User Weight	130kg (25st)
Product Weight	50kg



Tilted, safe seat day chairs

Ref. Singleton Mk1



Fabric choice

- Normal seat to back configuration with the whole body tilted i.e. the rear of the seat cushion is considerably lower than the front
- Fitted with four swivel, total lock brake castors for safety which must be braked when the user is sitting down/alighting or if the chair has been left unoccupied
- Wide pushing handle
- Wings fitted for snoozing
- Loose cushion with hygiene gap at rear for quick, effective cleaning
- Fitted with retractable heel support and supplied with a booster pad for those of smaller stature

Please note: These are heavy day chairs and are not for portering. They must only be used indoors on level and dry surfaces. They must not be moved when occupied.



Retractable heel support with booster pad

Height	1010mm (40")
Width	690mm (27")
Depth (Seat Retracted)	940mm (37")
Seat Height	580mm (23")
Seat Width	440mm (17")
Seat Depth	510mm (20")
Arm Height	690mm (27")
Back Height	660mm (26")
Max User Weight	160kg (25st)
Product Weight	40kg



Tilted, safe seat, semi-reclining day chairs

Ref. Singleton MkII



Fabric choice

- Available with fixed back or semi-reclining back with either standard or pressure relieving foam
- Normal seat to back configuration with the whole body tilted i.e. the rear of the seat cushion is considerably lower than the front
- Semi-reclining back can only be operated by the attendant
- Fitted with four swivel, total lock brake castors which must be braked when the user is sitting down or if the chair is unoccupied
- Wide pushing handle
- Spoon shaped buttoned backed for snoozing
- Retractable heel support and supplied with a booster pad for those of smaller stature

Please note: These are heavy day chairs and are not for portering. They must only be used indoors on level and dry surfaces. They must not be moved when occupied.



Semi-reclining with footrest and booster pad



Retractable heel support

Height	1010mm (40")
Width	690mm (27")
Depth (Seat Retracted)	940mm (37")
Seat Height	580mm (23")
Seat Width	440mm (17")
Seat Depth	510mm (20")
Arm Height	690mm (27")
Back Height	660mm (26")
Max User Weight	160kg (25st)
Product Weight	40kg



Tilted, safe seat day chairs

Ref. Fowey



Fabric choice

- Suitable for those with muscular problems and other special needs
- Gives maximum comfort and security
- Fully adjustable leg rest with self locking mechanism
- Four swivel, locking brake castors for safety
- Supplied complete with a removable, machine washable lambs wool cover which delivers pressure relieving capabilities
- A spare lambs wool cover may be purchased additionally if required



Height	970mm (38")
Width	640mm (25")
Length	1100mm (43")
Length with Leg Rest Extended	1500mm (59")
Max User Weight Guide	160kg (25st)
Product Weight	70kg